



## CyberScan Essentials

It is our company policy not to give any comments to the devices of our competitors. Below we list the features of our CyberScan system. It is left to the reader to compare and draw conclusions.

CyberScan is unique in several ways. As a genuine medical device (CE class IIa, TGA ARTG: DV-2012-DA-09896-3, FDA HCC D156965), CyberScan is ready for nearly everything a user may require.

The tough housing is made from aluminium rather than from cheap plastic. It can take a blow and is made for the harshest environment, like the Australian Outback.



The electronic parts inside the housing are developed and made in Germany; no cheap imported technique.

The colour touchscreen provides information about the status of the CyberScan and can be used to adapt the settings to ones needs (language, volume, connection etc.) and allows also access to the on-board operation system.

CyberScan is equipped with sufficient connection points to be ready for many accessories to come in the future without the need to change the design or hardware.





The software was developed by one of the finest software houses in Germany. Its use is very logical and intuitive and has a modular structure. That is the reason that this type of software surface is used by many other companies and perfectly adapted to their special needs.

We have chosen to use this software surface rather than developing a complete new one. But don't be fooled: everything in the background is one hundred per cent



CyberScan and includes (international edition)

- fully editable no-nonsense medical database
- already delivered with more than 107.000 informational signals in over 1.600 groups and sub-groups
- an absolutely unique algorithm to calculate the informational signals
- including the complete hierarchy of each signal in the calculations (therefore, one can differentiate between "normal" bacteria and antibiotic-resistant bacteria or determine whether an orange is an allergen or a beneficial food)


### Allergen




### Food




- templates for recurring tasks or for special tests (customised for a client or customised for specific tests like allergies or heavy metals etc.)



- creation of medical references – to create your own expert system



- biofeedback, bio resonance, radionics and colour- and sound tuning
- database for humans, animals (veterinarian use) and plants (agriculture)



- power supply from main power or battery pack
- high-speed software start (17 sec)
- easy exchange of groups and signals by built-in export and import tools
- on-line customer service for download of software and database updates
- on-line customer panel for exchange of experiences and information

If you compare this with other devices you'll know how to decide and – there should be only one choice.

AussiMed Pty. Ltd.  
[www.cyberscan.com.de](http://www.cyberscan.com.de)  
 Email: [mail@cyberscan.com.de](mailto:mail@cyberscan.com.de)  
 Phone: +49 160 977 14607